

ERISCROM

REFRIGERATED INCUBATOR FOR CHROMATOGRAPHS:

PRECISION AND CONTROL FOR YOUR LABORATORY

Multipurpose stainless steel refrigerated incubators with forced air circulation. With relative humidity control option. Ideal for applications with chromatographs, large stirrers, thermostatisation of instruments.



TEMPERATURE CONTROL

Standard temperature ranges are between +4°C and +55°C.



HUMIDITY CONTROL

Optional. Standard relative humidity range will be 30 to 90 %RH.



VOLUMES

Different volume options: 384, 718 or 1160 liters.



ECOLOGICAL GAS

Low GWP (value 3) R-290 gas in accordance with European F-GAS guidelines..



REINFORCED SHELF

Made of perforated stainless steel with fixed telescopic guides, and adjustable in height. Maximum load of 100 kg/m², specially designed for chromatographic applications.



CONTROLLER

Based on a microprocessor with 7" color touch screen, with electronic data logging system and graphical representation of the data. Indication of set and actual points.

